



www.bisonmedical.com



Celene

Multi CO₂ Laser





What is Cellene?

Cellene is the ultimate specialist in skin rejuvenation with both the CO₂ laser / RF-combined energy

It provides synergy effect by using the combined functions of CO₂ laser

As the name of Cell + Seline symbolizes, the multi RF & CO₂ energies shines the deteriorated cells/tissues back to the rejuvenation.



What consists to be all-in-one ?



What constitutes ?

CO2 Scanner

Unlike the conventional CO₂, it is adjustable in the 0.1mm Pitch. This is Fractional CO₂ laser that stimulates continuous collagen generation as well as helps to continuous skin rejuvenation by minutely destroying dermis and epidermis at small spot size.



Surgical scanning™

This unique Surgical Scanning™ can treat the previously untreatable bigger lesion smoothly by irradiating to max 10mm diameter. It can penetrate deeply and evenly. It can form the rectangle & oval shape, not only forming the typical square & circle shape.

Specifications

Laser system

Type	CO ₂ Laser
Wavelength	10,600nm, 650nm (Aiming Beam)
Max output power	60W
Delivery system	7-joint, fixed mirror articulated arm
Aiming beam	5mW, 650nm diode laser
Cooling system	Air
Electric requirements	100-230VAC, 50/60Hz, 1.5kVA
Dimension (W x D x H)	400 x 400 x 1880 mm
Weight	51 kg
Cut over current	6.3A over
Pulse width	0.2 - 350ms

Scan RF™ system

Scanning area	Up to 20 x 20 mm
RF type	Bipolar RF (Radiofrequency)
RF Frequency	2MHz
RF Power	Up to 50 Watt

RF type	
RF Frequency	2MHz
RF Power	Up to 50 Watt
Tips	Insulated and non-insulated tips

Hand pieces

Surgical mode	50, 100mm, Oral, Pharyngeal, Nasal Identical
Pixel mode	7 x 7, 9 x 9 pixels
CO ₂ mode	7 x 1 pixels

- **Rx only** Prescription only
- Indication for use, please refer attached from 3881